UAGB Nurse Triage AI Assistant - Quickstart Guide

Welcome!

The Nurse Triage Assistant helps you generate:

- Clinical Recommendations
- Patient Messages
- EMR Notes based on clinic protocols.

△ Step 1: Enter Patent Information

Fill out each of the 3 input boxes:

- Symptom Description (Example: "pink urine after stent removal")
- Relevant Recent Procedures or History (Example: "stent removal 2 days ago")
- Concerning Symptoms (Example: "mild dysuria, no fever, no clots"

% Step 2: Click Generate

Click the **Generate** button after filling all fields.

The system will produce:

- A Clinical Recommendation
- A Patient Message
- An EMR Note for charting

Step 3: Copy Outputs

- Click Copy All under each generated section.
- Paste the outputs into your patient chart, call message, or documentation system.

Important Reminders:

- Always **review** AI-generated notes before finalizing.
- Use clinical judgment to modify notes if needed.

• Escalate any uncertain situations to the supervising provider.

Quick Tips:

- Use simple, clinical language when filling fields.
- Focus on the most critical symptoms first.
- Save time by copying outputs directly into the EMR or patient portal messages.